

Thank you for considering Earthlite Lift tables. We offer choices on bases and tops, and we hope this step-by-step guide will help you selecting the right table for you. Please contact your local distributor for additional information.



1. Selecting base: select either Everest, Ellora, or Ellora Vista
 - Everest: Single pillar design provides stability, elegant look and more leg room (see Pic1 and Pic2)
 - Ellora: Unique scissor design offers wider height range (see Pic3 and Pic4)
 - Ellora Vista: Keeping wider height range, this unique curved legs and the hidden control box provides more sleek look, and this table offers the most working weight (see Pic5 and Pic6)

2. Selecting top
 - Flat top: No tilt (see Pic6 – headrest not included)
 - Tilt top: With back tilt – available in manual or power assist (see #3) (see Pic2 and Pic3)
 - Salon top: With back and leg tilt – available in manual or power assist (see #3)

3. Selecting adjustments (tilt and salon only)
 - Manual: Tilt adjustment is manual with 5 angles (20°, 30°, 50°, 55°, and 60°)
 - Power Assist: Tilt adjustment is non electric pneumatic with infinite adjustments between 0° and 56° by one touch lever

4. Selecting accessories (salon only)
 - Manual salon: Standard with salon pillow, neck roll, flex armrest, and footrest (see Pic4)
 - Power Assist: Option of with or without accessories (salon pillow, neck roll, and flex armrest) (see Pic1 and Pic5)

5. Selecting width: available in 28, 30, or 32 inches



Manual



Power Assist